



Recombinant Drosophila melanogaster Defense protein I (2)34Fc (I (2)34Fc)

Product Code	CSB-EP895094DLU
Abbreviation	I(2)34Fc
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V3Y3
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	YSD GAPKAACRDL TPQHGAKLQV TKPPYSISGP SHVRSQKLT LTLGGDEFLG FMIQARDGQN RVVGQFQVVD SVHSQTLDCS GKDDTITHLS AQKGKPLTGI TFDWIPPAGY KGNVKFMATV VQTGFVYWVG RVTKDIDVE
Source	E.coli
Target Names	I(2)34Fc
Protein Names	Recommended name: Defense protein I(2)34Fc
Expression Region	18-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.