



Recombinant *Drosophila melanogaster* MIP18 family protein CG30152 (CG30152)

Product Code	CSB-BP895108DLU
Abbreviation	CG30152
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V968
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLSYIKRKLS ESDSGVSSVA TVTSSCGGDS GRAGGTGSSE SGTGSSSASI SGRSQNADEL VRKTSQMSMD DEIAIFGEDA LLHELGYKNQ TELQETIYDL LRGIRDPEKP CTLEDLNVVY EDGIFVMPPT RSNVSVVRIE FNPTVPHCSL ATLIGLCIRV KVERGLPHNI KLDIYIKKGA HQTEEEINKQ INDKERIAAA MENPNLRDLV ENCIKDEE
Source	Baculovirus
Target Names	galla-1
Protein Names	Recommended name: MIP18 family protein CG30152
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.