



Recombinant Drosophila melanogaster Probable uridine-cytidine kinase (CG6364)

| | |
|--------------------------|---|
| Product Code | CSB-YP895117DLU |
| Abbreviation | CG6364 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9VC99 |
| Product Type | Recombinant Protein |
| Immunogen Species | Drosophila melanogaster (Fruit fly) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MQDLRNADTL RLPNGDGAAA NDEVKSPFLI GVAGGTASGK STVCKKIMEQ LGQAEMDHTQ RQVVSISQDS FYRELTPAEK AKAQKGLFNF DHPDAFNEEL MYSTLQNILK GHKVEIPSYD YRTNSLDFEN VLVYIPADV LFEGLVYFYF PKIRELFHMK LFDVTDSDTR LARRVPRDIN ERGRDLDAVL TQYMTFVKPA FEEFCSPKK FADVIPRGA DNTVAIDLIV HHIGEILATT NSAQHSNTVR VAASSMKRDH |
| Source | Yeast |
| Target Names | CG6364 |
| Protein Names | Recommended name: Probable uridine-cytidine kinase Short name= UCK EC= 2.7.1.48 Alternative name(s): Cytidine monophosphokinase Uridine monophosphokinase |
| Expression Region | 1-260 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |