



Recombinant *Drosophila melanogaster* Cysteine and histidine-rich domain-containing protein (CHORD)

Product Code	CSB-EP895118DLU
Abbreviation	CHORD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VCC0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEQCYNRGCG QLFDPQTNND ESCRHHPGEP FFHDAYKGWS CCNKKSVDFT EFLNIKGCTL AKHSNVKPPE PEKPKDES D KDEVIEVRAP IREALPRPPI DSPLTVIQPT VAPALKDMVF AVKTPAAQKS SDAIEVGTTC KNNGCTYSFT GNSSDFGECT YHPGVPIFHE GMKFWSCCQK RTSDFSQFMA QKGCTYGEHK WVKENDDKKV VQCRYDWHQT ATNVVMAIYA KKYDYSQSVI ELNPIRLHVN LVFPEQDNAR FDLLELRGI VNVSNASAHM YGTKVEIKLP KLEPGSWSNL NFPNKKLPVV KKSQVEEKKK QEESDEEFFD LDDIKAETSF RLSEMSMQSP NNLD
Source	E.coli
Target Names	CHORD
Protein Names	Recommended name: Cysteine and histidine-rich domain-containing protein Alternative name(s): Morgana
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.