



Recombinant *Drosophila melanogaster* Putative OPA3-like protein CG13602 (CG13602)

Product Code	CSB-YP895119DLU
Abbreviation	CG13602
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VCG2
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MVIGSFPVVGK IFIHGVKRIS KPFGDLLIWM GKHHPFIRRY IIIPPAQLYN TFEVRSKLRM LRLKQPRRIP PLSTPVATRL GADMLSEAFV FGIGAGLIYY ELSKCELLQM S
Source	Yeast
Target Names	CG43998
Protein Names	Recommended name: Putative OPA3-like protein CG13602
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.