



Recombinant Drosophila melanogaster Probable prefoldin subunit 3 (CG6719)

Product Code	CSB-BP895138DLU
Abbreviation	CG6719
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VGP6
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPKLPENQKT FAGIPEAVFL EEIDTFMSQP ENENCEKVLQ RLDEQHGKYR FMACNLEARR RKLKSQIPDL ERSLEMVNVL RKEDEERETQ FLLSDQVFIK TLVPPTKTVY LWLGASVMLE YPLDEAEALL NQNITSAVGN LKSVEHDQDF LRDQITTEV NMARVYNWGV KKRQAATKTT ATTPS
Source	Baculovirus
Target Names	mgr
Protein Names	Recommended name: Probable prefoldin subunit 3
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.