



# Recombinant *Drosophila melanogaster* JmjC domain-containing protein 4 homolog (CG7200)

<b>Product Code</b>	CSB-EP895150DLU-B
<b>Abbreviation</b>	CG7200
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9VJ97
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASEERDGDV LLQLHPPEEVV ETDPLPEEI ERCSGLDYND FFWRYMHKNI PVIIANVSND WECQNWTVGQ SSPEERDLNS NPSASSINFD YLTKISDGP VPVANCNSSY FNSHTKLELN FHDYLAKWRS SIESQSSAAW TSAEVNSNVA PASGDNLYLK DWHLAAQMPG YNFYKVPKYF ASDWLNEQLI QQGKDDYRFV YMGPKNSWTS YHADVFGSFS WSTNIVGLKK WLIMPPGEEL KLNDRLGNLP FSIDEKMLDE HNVRYTINQ RANEAVFVPS GWFHQVWNLDT DTISVNHWNF NGCNISMVWQ NLKNNLKAVC NEISDCQQMD NFEAHCQTML RASFGINYLD FIELLEFIAA RRLAEGTVAT KFLLFDSYTM NDYHVQYDLE CLWKITRTL EDPTIQCSPL QLEDRCQGLL ARLEF
<b>Source</b>	E.coli
<b>Target Names</b>	CG7200
<b>Protein Names</b>	Recommended name: JmjC domain-containing protein 4 homolog Alternative name(s): Jumonji domain-containing protein 4 homolog
<b>Expression Region</b>	1-425
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.