



Recombinant Drosophila melanogaster Probable elongation factor 1-delta (eEF1delta)

Product Code	CSB-EP895157DLU
Abbreviation	eEF1delta
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VL18
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MKVEALDKFW ADKSRDYLAE KRFYEGPQKV TDRSHYSPLV SEIAKAREHI QNSLEKIDGV TLDDGLNSEL AKRLAQLEGE HKELKTQVSL LNELLTATVK RLETQLKLTN GVSKEPEVEA KKPEANDDDD DVDLFGSDSE EEDGEAARIR EERLAAYAAK KAKKVQIIAK SNIILDVKPW DDETDLKVME TEIRKITQDG LLWGASKFVP VAFGIQKLSI SCVVEDDKVS IDWLTEEIEK LEDFVQSVDI AAFNKI
Source	E.coli
Target Names	eEF1delta
Protein Names	Recommended name: Probable elongation factor 1-delta Short name= EF-1-delta
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.