



Recombinant *Drosophila melanogaster* SOSS complex subunit B homolog (CG5181)

Product Code	CSB-BP895163DLU
Abbreviation	CG5181
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VM17
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MYNVECIPIK DIKPGLKNIN VIFIVLEVGV ATVTKENREV RNFKVGDPTA CINVSIWDEP GKLIAPGDIV RLTKGYASIW RHCLTLYSGK NGEVFKIGEY CMVFNESVNM SEPKRAEQQA VANPAATPAG LPAGGGAPGL PAKGGATGIP QPAVAAAPGA PATQSAVTTA PAAAPAIAPQ TTTKPGTRGG RGGGGRGGLK GERR
Source	Baculovirus
Target Names	CG5181
Protein Names	Recommended name: SOSS complex subunit B homolog
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.