



Recombinant *Drosophila melanogaster* Golgi to ER traffic protein 4 homolog (CG9853)

Product Code	CSB-BP895169DLU
Abbreviation	CG9853
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VN19
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSAAETGPPIA SGSAVSTGQR GVSRLAKLS QSLAGGEFYE AHMMYRTLYF RYTAQKRYQD CLELLFDGAQ QLIAKEQESS AADLCLLLVD TLEKRGPAE DTDNFLWVPR LGALIRGLNS ATVERETLIQ RTIKWSTALH GQYGHPVLHK LIAHVFWTEG NIESARHHYL LCQDGS LCGR VLIEISQSRG FQGEMDLFLV QAVLQQLSLK DRKTAEDTFT EYTRYHAKLL KHEFPYKEPL VNFLYFLFRL IDAKRVAGFR ALRKLYDPSL KRDT SFLKYV AKIGVLYFDE HPEAAHSGPP GLGGMFGDIF NRLMAGFEED DADDHDTTLR RQENNELE
Source	Baculovirus
Target Names	CG9853
Protein Names	Recommended name: Golgi to ER traffic protein 4 homolog
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.