



Recombinant Drosophila melanogaster Eukaryotic translation initiation factor 3 subunit F-1 (eIF3-S5-1)

Product Code	CSB-YP895170DLU
Abbreviation	eIF3-S5-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VN50
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSALNLTVRV HPVVLQVVD AFERRNADSH RVIGTLLGSV DKGVVEVTNC FCVPHKEHDD QVEAELSYAL DMYDLNRKVN SNESVVGWWA TGNDVTNHSS VIHEYAREC NNPVHLTVDT SLQGGRMGLR AYVCIQLGVP GGKSGCMFTP IPVELTSYEP ETFGLKLLQK TVGVSPAHRP KTVPPMLDLA QISEASTKLQ SLLDLILKYV DDVIAHKVTP DNAVGRQLLD LIHSVPHMTH EQFTQMFNAN VRNLLLVTIL SQLIKTQLQL NEKLTFLPTA
Source	Yeast
Target Names	eIF3f1
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit F-1 Short name= eIF3f-1 Alternative name(s): Eukaryotic translation initiation factor 3 subunit 5-1
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.