



Recombinant *Drosophila melanogaster* KRR1 small subunit processome component homolog (dbe)

Product Code	CSB-YP895176DLU
Abbreviation	dbe
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VPU8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSESEAEETK ISTEPVDNAW SMKIPAFRQE DNPFGMVEES SFATLFPKYR ERYLKEVWPL VEQCLAEHHL KAELDLMEGS MVVKTSRKTW DPYIIKARD MIKLMARSVP FEQAKRVLQD DIGCDIIKIG NLVHKKEKFV KRRQRLIGPN GATLKSIELL TDCYVLVQGN TVSALGPYKG LQQVRDIVLE TMNNVHPIYN IKALMIKREL MKDPRLANED WSRFLPKFKN KNISKRKQPK VKKQKKEYTP FPPSQPESKV DKQLASGEYF LNQEQQQAKR NQERTEKQKE AAKRQDERRN KDFVPPTEES AASSRKKEDG SSSSKVDVKA LKAKLIKANK KARSS
Source	Yeast
Target Names	dbe
Protein Names	Recommended name: KRR1 small subunit processome component homolog Alternative name(s): KRR-R motif-containing protein 1 Protein dribble
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.