



Recombinant *Drosophila melanogaster* THO complex subunit 6 (thoc6)

Product Code	CSB-MP895192DLU
Abbreviation	thoc6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VTV1
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MKQKDLKRAY NNVLAQAISG SKQYLFAGNL FGDIFVLRİK ELDKGSEPP GKLKIFPQGS DVDINYLA FH RDFLIVGAVG LIYGLEWNEE EESLATKRSW EVKIPMQVDA VEPDVNSMW LDSENSILFA GCGDGVİYQV SLEDGRIQRE YRGHTDYVHS VVG NANGQIF SGAEDGTVRV WSTKQQQHTS MLEPYKNPNL LRPDWGKWIG AVAVNEDWLL CGGGPKASIF HLRSMESTCV FSFPGRVHLC DFVDDCVLIG GEHNHVQSYT LNGVLQANIP VEHTACYS AV WQTSPİKFIS IAGFSNKLHI LKDFRFLDSK İEMYGNVEGE ENDPEEEDEL IEKPLTVEAA
Source	Mammalian cell
Target Names	thoc6
Protein Names	Recommended name: THO complex subunit 6
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.