



# Recombinant Drosophila melanogaster Probable cytochrome P450 308a1 (Cyp308a1)

<b>Product Code</b>	CSB-MP895206DLU
<b>Abbreviation</b>	Cyp308a1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9VWR2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila melanogaster (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLPLVLFILL AATLLFWKWQ GNHWRRLLGLE APFGWPLVGN MLDFALGRRS YGEIYQEIYT RNPGLKYVGF YRLFNEPAIL VRDQELLRQI LVGRNFADCA DNAVYVDHQR DVLASHNPFI ANGDRWRVLR ADLVPLFTPS RVRQTLPHVA RACQLLRDQV PLGRFEAKDL ATRYTLQVVA SAIFGLDAHC LGIHMVAHE PSRWLEWLAP LFQPSVWVSL ETMSLLHTPR LGRLIGHRYV PLPLQHWFRE LVEARSGGDN LLQWLAESKR GLGKEELAGH ATTLLEGYE TSAMLLAFAL YELALNEDAQ RRLHIELDEV AQRHAGNLID PVALGELRYS EAALLEALRL HPAMQALQKR CTKTFTLPDQ KSGASSELKV HLGTVLVLPV QAIHLDPALY PAPNQFRPER FLNQPPMGCR FLGFGAGPRM CPGMRLGLLQ TKAALTLLQ DHCVQLADED QCRVEVSPLT FLTASRNGIW LSKRRTRRY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Cyp308a1
<b>Protein Names</b>	Recommended name: Probable cytochrome P450 308a1 EC= 1.14.-.- Alternative name(s): CYPCCCVIIIA1
<b>Expression Region</b>	1-490
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.