



Recombinant *Drosophila melanogaster* 60S ribosomal protein L5 (RpL5)

Product Code	CSB-MP895242DLU
Abbreviation	RpL5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W5R8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGFVKVVKNK QYFKRYQVKF RRRREGKTDY YARKRLTFQD KNKYNTPKYR LIVRLSNKDI TVQIAYARIE GDRVVCAAYS HELPKYGIQV GLTNYAAAYC TGLLVARRVL NKLGLDSLVA GCTEVTGEEF NVEPVDDGPG AFRCFLDVGL ARTTTGARVF GAMKGAVDGG LNIPHSVKRF PGYSAETKSF NADVHRAHIF GQHVADYMRS LEEDEESFK RQFSRYIKLG IRADDLEDIY KKAHQAIRND PTHKVTAKKS SAVTKKRWNA KKL TNEQRKT KIAAHKAAYV AKLQSETEA
Source	Mammalian cell
Target Names	RpL5
Protein Names	Recommended name: 60S ribosomal protein L5
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.