



Recombinant Mouse STIP1 homology and U box-containing protein 1 (Stub1)

Product Code	CSB-MP895285MO
Abbreviation	Stub1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WUD1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MKGKEEKEGG ARLGTGGGGS PDKSPSAQEL KEQGNRLFVG RKYPEAAACY GRAITRNPLV AVYYTNRALC YLKMQQPEQA LADCRRALEL DGQSVKAHFF LGQCQLEMES YDEAIANLQR AYSLAKEQRL NFGDDIPSAL RIAKKKRWNS IEERRIHQES ELHSYLTRLI AAERERELEE CQRNHEGHED DGHIRAQQAC IEAKHDKYMA DMDELFSQVD EKRRKRDPD YLCGKISFEL MREPCITPSG ITYDRKDIEE HLQRVGHFDP VTRSPLTQEQ LIPNLAMKEV IDAFISENGW VEDY
Source	Mammalian cell
Target Names	Stub1
Protein Names	Recommended name: STIP1 homology and U box-containing protein 1 EC= 6.3.2.- Alternative name(s): Carboxy terminus of Hsp70-interacting protein E3 ubiquitin-protein ligase CHIP
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.