



Recombinant Rat Insulin-like 3 (Insl3), partial

Product Code	CSB-MP895287RA
Abbreviation	Insl3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WUK0
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Source	Mammalian cell
Target Names	Insl3
Protein Names	Recommended name: Insulin-like 3 Alternative name(s): Leydig insulin-like peptide Short name= Ley-I-L Relaxin-like factor Cleaved into the following 2 chains: 1. Insulin-like 3 B chain 2. Insulin-like 3 A chain
Expression Region	-
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This protein is an insulin like hormone produced mainly in gonadal tissues in males and females. Studies of the mouse counterpart suggest that this gene may be involved in the development of urogenital tract and female fertility. It may also act as a hormone to regulate growth and differentiation of gubernaculum, and thus mediating intra-abdominal testicular descent. The mutations in this gene may lead to, but not a frequent cause of, cryptorchidism.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.