



# Recombinant Rat Segment polarity protein dishevelled homolog DVL-1 (Dvl1)

<b>Product Code</b>	CSB-EP895301RA
<b>Abbreviation</b>	Dvl1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9WVB9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAETKIIYHM DEEETPYLVK LPVAPERVTL ADFKNVLSNR PVHAYKFFFK SMDQDFGVVK EEIFDDNAKL PCFNGRVVSWS LVLAEGAHS AGSQGTDSHT DLPPPLERTG GIGDSRPPSF HPNVASSRDG MDNETGTESM VSHRRERARR RNRDEAARTN GHPRGDRRRE LGLPPDSAST VLSSELESSS FIDSDEEDNT SRLSSSTEQS TSSRLIRKHK CRRRKQRLRQ TDRASSFSSI TDSTMSLNII TVTLNMERHH FLGISIVGQS NDRGDGGIYI GSIMKGGAVA ADGRIE PGDM LLQVNDVNFE NMSNDDAVRV LREIVSQTGP ISLTVAKCWD PTPRSYFTIP RADPVRPIDP AAWLSHTAAL TGALPRYGTS PCSSAITRTS SSSLTSSVPG APQLEEAPLT VKSDMSAIVR VMQLPDSGLE IRDRMWLKIT IANAVIGADV VDWLYTHVEG FKERRARKY ASSMLKHGFL RHTVNKITFS EQCYVFGDL CSNLASLNLN SGSSGASDQD TLAPLPHPSV PWPLGQGYPY QYPGPPPCFP PAYQDPGFSY GSGSAGSQQS EGSKSSGSTR SSHRTPGREE RRATGAGGSG SESDHTVPSG SGSTGWWERP VSQLSRGSSP RSQASAVAPG LPPLHPLTKA YAVVGGPPGG PPVRELA AVP PELTGSRQSF QKAMGNPCEF FVDIM
<b>Source</b>	E.coli
<b>Target Names</b>	Dvl1
<b>Protein Names</b>	Recommended name: Segment polarity protein dishevelled homolog DVL-1 Short name= Dishevelled-1 Alternative name(s): DSH homolog 1
<b>Expression Region</b>	1-695
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	DVL1, the human homolog of the Drosophila dishevelled gene (dsh) encodes a cytoplasmic phosphoprotein that regulates cell proliferation, acting as a transducer molecule for developmental processes, including segmentation and



neuroblast specification. DVL1 is a candidate gene for neuroblastomatous transformation. The Schwartz-Jampel syndrome and Charcot-Marie-Tooth disease type 2A have been mapped to the same region as DVL1. The phenotypes of these diseases may be consistent with defects which might be expected from aberrant expression of a DVL gene during development.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.