



# Recombinant Mouse EH domain-containing protein 1 (Ehd1)

<b>Product Code</b>	CSB-EP895306MO-B
<b>Abbreviation</b>	Ehd1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9WVK4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFSWVSKDAR RKKEPELFQT VAEGLRQLYA QKLLPLEEHY RFHEFHSPAL EDADFDNKPM VLLVGQYSTG KTFIRHLIE QDFPGMRIGP EPTTDSFIIV MHGPTGVVP GNALVVDPRR PFRKLNAFGN AFLNRFMCAQ LPNPVLDSIS IIDTPGILSG EKQRISRGYD FAAVLEWFAE RVDRIILLFD AHKLDISDEF SEVIKALKNH EDKIRVVLNK ADQIETQQLM RVYGALMWSL GKIINTPEVV RVYIGSFWSH PLLIPDNRKL FEAEEQDLFK DIQSLPRNAA LRKLNLIKR ARLAKVHAYI ISSLKKEMPN VFGKESKKKE LVNNLGEIYQ KIEREHQISS GDFPSLRKMQ ELLQTQDFSK FQALKPKLLD TVDDMLANDI ARLMVMVRQE ESLMPSQAVK GGAFDGTMNG PFGHGYGEGA GEGIDDVEWV VGKDKPTYDE IFYTLSPVNG KITGANAKKE MVKSKLPNTV LGKIWKLADV DKDGLLDDEE FALANHLIKV KLEGHELPAD LPPHLIPPSK RRHE
<b>Source</b>	E.coli
<b>Target Names</b>	Ehd1
<b>Protein Names</b>	Recommended name: EH domain-containing protein 1 Alternative name(s): PAST homolog 1 Short name= mPAST1
<b>Expression Region</b>	1-534
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene belongs to a highly conserved gene family encoding EPS15 homology (EH) domain-containing proteins. The protein-binding EH domain was first noted in EPS15, a substrate for the epidermal growth factor receptor. The EH domain has been shown to be an important motif in proteins involved in protein-protein interactions and in intracellular sorting. This protein is thought to play a role in the endocytosis of IGF1 receptors.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.