



Recombinant Trehalose-phosphatase (gob-1)

Product Code	CSB-BP895482CXY
Abbreviation	gob-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XTQ5
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MFTFRLSTDL SLKLAKFPPG RAGSRPDCPM NCEKESQMTI ASQSIEDFKE CLYQMQEARK SVTNEILETG HIKADQVQIF KSTLEEMNDE RTSKNHIRDI HSRGTTFGIN IQDEIKGLQK DHHFLDAFAV ESDKENNSFA NVLKLCDLPG LLSKFVDDEI RFEKEVAECK AFLMDLIDTS TTGGIKPLFI TDWDGTMKDY CSQYATNLQP AYSIVMGVF SRLFTRAFV LTAGPLRHPG ILDLTALPIN GPVLFSGSWG REWWLGGRRI VHDDGIPEEG SVAIGQLCEQ LDEILHEGEG VQFALVGSGV QRKVDRLTLG VQTVFKQVPE DLSARYIDAV RERIHRVDPN SQYLVLENC S PLEIEVCVHS SGAVWNKGDG VAALVESLHD SLKVGKVCVA GDTASDV PML KKAADENPEN VRALFVNINK QLQENITNIV GDAKRVCFIS SPDVAHAAFA QIIEFSG
Source	Baculovirus
Target Names	gob-1
Protein Names	Recommended name: Trehalose-phosphatase EC= 3.1.3.12 Alternative name(s): Trehalose-6-phosphate phosphatase Short name= TPP
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.