



Recombinant Drosophila melanogaster Probable cytochrome P450 6a13 (Cyp6a13)

Product Code	CSB-MP895568DLU
Abbreviation	Cyp6a13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V4U9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLTLLVLVFT VGLLLYVKLR WHYSYWSRRG VAGERPVYFR GNMSGLGRDL HWTDLNLRiy RKFRGVERYC GYFTFMTKSL FIMDLELIRD IMIRDFSSFA DRGLFHNVRD DPLTGNLLFL DGPEWRWLRQ NLTQVFTSGK MKFMFPNMVE VGEKLTQACR LQVGEIEAKD LCARFTTDVI GSCAFGLECN SLQDPESQFR RMGRSVTQEP LHSVLVQAFM FAQPELARKL RFRLFRPEVS EFFLDTVRQT LDYRRRENIH RNDLIQLLME LGEEGVKDAL SFEQIAAQUAL VFFLAGFDTS STTMSFCLYE LALNPDVQER LRVEVLAVLK RNNQKLTYS VQEMPYLDQV VAETLRKYPI LPHLLRRSTK EYQIPNSNLI LEPGSKIIP VHSIHDPPEL YPDPEKFDPS RFEPEEIKAR HPFAYLPFGE GPRNCIGERF GKLQVKVGLV YLLRDFKFSR SEKTQIPLKF SSRNFLISTQ EGVHLRMEGL ERP
Source	Mammalian cell
Target Names	Cyp6a13
Protein Names	Recommended name: Probable cytochrome P450 6a13 EC= 1.14.-.- Alternative name(s): CYPVIA13
Expression Region	1-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.