



Recombinant *Drosophila melanogaster* Guanosine-3',5'-bis (diphosphate) 3'- pyrophosphohydrolase MESH1 (Mesh1)

| | |
|--------------------------|---|
| Product Code | CSB-MP895585DLU |
| Abbreviation | Mesh1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9VAM9 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila melanogaster</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MATYPSAKFM ECLQYAAFKH RQRRKDPQE TPYVNHVINV STILSVEACI TDEGVLMAAL LHDVVEDTDA SFEDVEKLFV PDVCGLVREV TDDKSLEKQE RKRLQIENAA KSSCRAKLIK LADKLDNLRD LQVNTPTGWT QERRDQYFVW AKKVVDNLRG TNANLELKLD EIFRQRGLL |
| Source | Mammalian cell |
| Target Names | Mesh1 |
| Protein Names | Recommended name: Guanosine-3',5'-bis(diphosphate) 3'- pyrophosphohydrolase MESH1 EC= 3.1.7.2 Alternative name(s): Metazoan SpoT homolog 1 Short name= Mesh1 Penta-phosphate guanosine-3'- pyrophosphohydrolase Short name= (ppG) |
| Expression Region | 1-179 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |