



# Recombinant Drosophila melanogaster E3 ubiquitin-protein ligase RING1 (Sce)

<b>Product Code</b>	CSB-MP895587DLU
<b>Abbreviation</b>	Sce
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9VB08
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila melanogaster (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSLDPAPNK TWELSLYELQ RKPQEVITDS TEIAVSPRSL HSELMCPICL DMLKKTMTTK ECLHRFCSDC IVTALRSGNK ECPTCRKKLV SKRSLRADPN FDLLISKIYP SREEYEAIQE KVMAKFNQTQ SQQALVNSIN EGIKLQSQNR PQRFRTKGGG GGGGGGGNGN GAANVAAPPA PGAPTAVGRN ASNQMHVHDT ASNDSNSNTN SIDRENDRPG HSGTSAASAI TSASNAAPSS SANSASTSA TRMQVDDASN PPSVRSTPSP VPSNSSSSSKP KRAMSVLTSE RSEESDSQ MDCRTEGDSN IDTEGEGNGE LGINDEIELV FKPHPTMSA DNQLIRALKE NCVRYIKTTA NATVDHLSKY LAMRLTLDLG ADLPEACRVL NFCIYVAPQP QQLVILNGNQ TLHQVNDKFW KVNKPMEMYY SWKKT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Sce
<b>Protein Names</b>	Recommended name: E3 ubiquitin-protein ligase RING1 EC= 6.3.2.- Alternative name(s): Sex comb extra protein dRING protein dRING1
<b>Expression Region</b>	1-435
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.