



Recombinant *Drosophila melanogaster* Uncharacterized protein CG5902 (CG5902)

Product Code	CSB-BP895592DLU
Abbreviation	CG5902
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VCF0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MSNSHYNNYQ QQQPHSSNGD PEYQHQQMVH QPQRFSNGHG MKTVAVPDMC LFCFEVLDCE LNNVDGPSVP VFSNDAYPLF VTKIGRDKR LRGCIGTFSA MELHHGLREY ALTSAFKDSR FAPISRDELP RLTVSVSILQ NFEEAQGHLD WQLGVHGIRI EFLTERGCKR TATYLPQVAT EQGWDQLQTI DSELLRKGGYR AAITPETRKS IKLTRYRSQE IQMHYKEYRE YQERRAQFGK VQC
Source	Baculovirus
Target Names	CG5902
Protein Names	Recommended name: Uncharacterized protein CG5902
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.