



Recombinant *Drosophila melanogaster* UPF0585 protein CG18661 (CG18661)

Product Code	CSB-EP895628DLU
Abbreviation	CG18661
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VLF6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MSLTKRLLHI HPRCSFSKRS HPSADRNSQP ISEALLSQVD KTTENLQLE IASGSGQHAG FLAPLLPNIT FQPTEYERNQ FGSIAAYAGD CPTRNIRPPF HVDITRDPGD WEAPPAPASY DYMFNSNMMH ISPWSCSIGL FRAAGQLLKK GGRMFTYGPY AQDGILVPQS NVDFDRSLRQ RDASWGVVDI KDLKVLAEN GLQLEKLVKM PSNNKFLTWL KL
Source	E.coli
Target Names	CG18661
Protein Names	Recommended name: UPF0585 protein CG18661
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.