



Recombinant *Drosophila melanogaster* Ras-related protein Rab-26 (Rab26)

Product Code	CSB-BP895640DLU
Abbreviation	Rab26
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VP48
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MASTAVGLGG GEGDPGAGGP PAGSAHPDDA SSMSDDVFED AETTQARIEE LRRRPFGDGS YNPPAAPASV SASITTTTTQ QQQQHNNPSH HHQSSHHQPS HHHHHHHHSQ LSLTGSHHYH DDAIMAPVQR SATGYPGYRP SREAMQMYAY GTDDYDDDDYN DGWRSYRYDE VDMHPAPSNA HQQPFDDTVN HKTILLGDSG VGKTSFLVKY NTGEFRLGSF SATVGIALTN KVVVVDGTRV KLQIWDTAGQ ERFRSVTHAY YRDAHALLL YDVTNKTTYD NIRAWLGEIR EYAQEDVVIV LIGNKADCSG SERQVKREDG ERLGREHNVP FMETSAKTGL NVELSFTAVA RQLKSRGYEH GDDGKFNVDH FVRDNTKARS VCAQC
Source	Baculovirus
Target Names	Rab26
Protein Names	Recommended name: Ras-related protein Rab-26
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.