



Recombinant *Drosophila melanogaster* Serine protease persephone (psh)

Product Code	CSB-MP895671DLU
Abbreviation	psh
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VWU1
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	AVTVGRAC KVTDTMPGIC RTSSDCEPLI DGYIKSGVLT LNDVPSCGLG AWGEIFCCPT KPCCDNSTIT SVSTSSTTST KAPMTSGRVD VPTFGSGDRP AVAACKKIRE RKQQRSGNQL VIHIVGGYPV DPGVYPHMAA IGYITFGTDF RCGGSLIASR FVLTAAHCVN TDANTPAFVR LGAVNIENPD HSYQDIVIRS VKIHPQYVGN KYNDIAILEL ERDVVETDNI RPACLHTDAT DPPSNSKFFV AGWGVLVNVT RARSKILLRA GLELVPLDQC NISYAEQPGS IRLKQGVID SLLCAIDQKL IADACKGDSG GPLIHELNVE DGMYTIMGVI SSGFGCATVT PGLYTRVSSY LDFIEGIVWP DNRV
Source	Mammalian cell
Target Names	psh
Protein Names	Recommended name: Serine protease persephone EC= 3.4.21.-
Expression Region	23-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.