



Recombinant *Drosophila melanogaster* Protein NipSnap (Nipsnap)

Product Code	CSB-YP895674DLU
Abbreviation	Nipsnap
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VXK0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLKLRNLLAV GKSNNNAVRS LSTTPSRNDS ESWFSKLLVR KIEPTKESHS RMLSDKEIY ALHTHNVRPD SMGSYLNNYK TTVALINEKK ANLSCELVAS WTVQVGDMQ CLHLWKYTGK FEKIDQAKED LWNDPEYLSL MQERSKFLRS RHLQYLLAFS YWPQIASRTG KNIYEMRSYR LTPGTMIEWG NNWARAINYR KHNNEAFAGF FSQIGRLYNV HHIWCYKSLQ DRKETREAAW RSPGWDECVA YTVPLIREMH CRVLAPTEFS PSQ
Source	Yeast
Target Names	Nipsnap
Protein Names	Recommended name: Protein NipSnap
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.