



Recombinant *Drosophila melanogaster* Neuropeptide-like 1 (Nplp1)

Product Code	CSB-EP895687DLU
Abbreviation	Nplp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W0W6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	SVAALAAQG LLNAP
Source	E.coli
Target Names	Nplp1
Protein Names	Recommended name: Neuropeptide-like 1 Cleaved into the following 9 chains: 1. MTYamide peptide 2. IPNamide peptide 3. NAP peptide 4. NPLP1-1 5. NPLP1-2 6. NPLP1-3 7. NPLP1-4 8. NPLP1-5 9. NPLP1-6
Expression Region	102-115
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.