



Recombinant *Drosophila melanogaster* UPF0670 protein CG4666 (CG4666)

Product Code	CSB-EP895698DLU-B
Abbreviation	CG4666
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W440
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	FI RCVFTVFAGR LFQRKRKVTDTTTIYGLCTS QDVDIFIRHM NNARYLRELD FARFHFYALT GLYERIRDRR GGAVQGASSV RYRRTIPIFH PYKIQTCLIW WDDKAIYLEQ QFITLSDGFV RAVAMSKQNI TNCNVLEVLK TYPETAQRPE KPEELKLWLD AIELSSQKLR KDK
Source	E.coli
Target Names	CG4666
Protein Names	Recommended name: UPF0670 protein CG4666
Expression Region	19-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.