



Recombinant *Drosophila melanogaster* Sel1 repeat-containing protein 1 homolog (CG13865)

Product Code	CSB-EP895704DLU-B
Abbreviation	CG13865
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W5N0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAYDLKKESD VKEYVEKLGV EYRFGCYSEK KPEACHLLGD YLEGIKKDFE KASKVYKSTC DDYGYAKSCY KYGNYSFLGK GKSGSKGNPQ VAYEYYEKGC NLNDSACLH SGLLLVSKSM PREIDWNVPK GLEFLTKSCD LNNATACFYL SGMHISGVQK KADQSAVTAS SGSGTSSPPA GQPPLKDSY IVLKDMKKAF QFAHKACELR NMYACANLSQ MYARGDGIK NEKEAEKYKK LALEMQDEVK KQHDTLGFQQ GVGMPN
Source	E.coli
Target Names	CG13865
Protein Names	Recommended name: Sel1 repeat-containing protein 1 homolog Alternative name(s): Beta-lactamase hcp-like protein CG13865
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.