



Recombinant Rat 60S ribosome subunit biogenesis protein NIP7 homolog (Nip7)

Product Code	CSB-BP895759RA
Abbreviation	Nip7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WV50
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MRPLTEEETR VMFEKIAKYI GENLQLLVDR PDGTYCFRLH NDRVYYVSEM ILKLAANISG DKLVSLGTCF GKFTKTHKFR LHVTALDFLA PYAKYKWWIK PGAEQSFLYG NHVLKSGLGR ITENTSQYQG VVVYSMADVP LGFGVAAKST QDCRKVDPMA IVVFHQADIG EYIRHEETLT
Source	Baculovirus
Target Names	Nip7
Protein Names	Recommended name: 60S ribosome subunit biogenesis protein NIP7 homolog Alternative name(s): PEachy kDa93
Expression Region	1-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.