



Recombinant Arabidopsis thaliana Ras-related protein RABG3e (RABG3E)

Product Code	CSB-EP895906DOA
Abbreviation	RABG3E
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XI98
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MPSRRRTLLK VIILGDSGVG KTSLMNQYVN KKFSNQYKAT IGADFLTKEV QFEDRLFTLQ IWDTAGQERF QSLGVAFYRG ADCCVLVYDV NSAKSFEDLN NWREEFLIQA SPSPENFPF VVIGNKIDVD GGSSRVVSEK KARAWCASKG NIPYYETSAK VGTNVEDAFL CITTNAMKSG EEEEMYLPDT IDVGTSNPQR STGCEC
Source	E.coli
Target Names	RABG3E
Protein Names	Recommended name: Ras-related protein RABG3e Short name= AtRABG3e Alternative name(s): Ras-related protein Rab74 Short name= AtRab74
Expression Region	1-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.