



Recombinant Bovine Alpha-enolase (ENO1)

Product Code	CSB-EP895935BO-B
Abbreviation	ENO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XSJ4
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	SILKVHARE IFDSRGNPTV EVDLFTAKGL FRAAVPSGAS TGIYEALRLR DNDKTRYMGK GVS KAVEHIN KTIAPALVSK KLVNVEQEKI DKLMIEMDGT ENKSKFGANA ILGVSLAVCK AGAVEKGVPL YRHIADLAGN AEVILPVPAP NVINGGSHAG NKLAMQEFMI LPVGAENFRE AMRIGAEVYH NLKNVIKEY GKDATNVGDE GGFAPNILEN KEALELLKNA IGKAGYSDKV VIGMDVAASE FYRSGKYDLD FKSPDDPSRY ITPDELANLY KSFIRDYPVV SIEDPFDQDD WEAWQKFTAS AGIQVVGDDL TVTNPKRIAK AVSEKSCNCL LLKVNQIGSV TESLQACKLA QSNWGWVMVS HRSGETEDTF IADLVVGLCT GQIKTVAPCR SERLAKYNQI LRIEEELGSK AKFAGRSFRN PLAK
Source	E.coli
Target Names	ENO1
Protein Names	Recommended name: Alpha-enolase EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase Enolase 1 HAP47 Non-neural enolase Short name= NNE Phosphopyruvate hydratase
Expression Region	2-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes one of three enolase isoenzymes found in mammals; it encodes alpha-enolase, a homodimeric soluble enzyme, and also encodes a shorter monomeric structural lens protein, tau-crystallin. The two proteins are made from the same message. The full length protein, the isoenzyme, is found in the cytoplasm. The shorter protein is produced from an alternative translation start, is localized to the nucleus, and has been found to bind to an element in the c-myc promoter. A pseudogene has been identified that is located on the other arm of the same chromosome.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.