



Recombinant Probable phosphomannomutase (F52B11.2)

Product Code	CSB-EP895948CXY-B
Abbreviation	F52B11.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XUE6
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTSPASSRTI LVFDVDGTLT AARQTITPEM REFLIEARKR VPLAIVGGSD FKKITEQLAD HDKDLLLSLF DYTFSENGLY GFKGTEPYPV QSIQKAIGDA KLQELINFAL RYMSDIQLPV KRGNFVEFRN GMINLSPIGR SCSQEERMQF VEFDKKGIR QKFTEQLREK FGQYGLQFAI GGQISVDVFP TGWDKTFCLQ YLVPDFDTIH FFGDKTAPGG NDHEIFADER TVGHTVEGPE DTRKHVENVL KELD
Source	E.coli
Target Names	F52B11.2
Protein Names	Recommended name: Probable phosphomannomutase Short name= PMM EC= 5.4.2.8
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.