



Recombinant Mediator of RNA polymerase II transcription subunit 29 (mdt-29)

Product Code	CSB-BP895949CXY
Abbreviation	mdt-29
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XUS2
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSGQGPPPSM TPQQQHMLMQ QQQQQQMMRQ QQMQQQQMQQ QRQQQIQQA QQPYQRARTP QMVPQGGSPG GAHLQMHPHL QPQQQMQRPS PLVGAQLQAP SSVPTAANPT TPQMMQQQMG MNQPMSLPPT HVSRRPGSVAP PSSVPVNLQH TSGAPGLPGS QMEHQYPMHL QPQQQTLSRP GSQQSQHIQQ PGSIQRPGSV LAPGSIPPGG PASQTGPQSI QVFGPGSVQP PGSTQAPSSV QPASTFNPGS IQAPASQQPP ASVQPPPSAA SGSTVAGAQS SKEPLKPNEE QIRMVQDPVD LVRNLVQKDL RNSLVEMNKR GAELVRQKQE GEVNEDGKAQ YSRATNDFHA VCDEIDRTLTL TVLETAKQLS KLDKVFFDRS SRDLDEVMV NSVQNFVDNT EIVQRMFDDT IGSVTSSMES MRRRQKKWED QHKNDDVEMM E
Source	Baculovirus
Target Names	mdt-29
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 29 Alternative name(s): Mediator complex subunit 29
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.