



Recombinant Protein smg-9 (smg-9)

Product Code	CSB-EP895955CXY
Abbreviation	smg-9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9XWJ1
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MKKVEILKTS RPSSAGGAAR PSTASPTHGA PKIAIKTRPV ADDVAPTAAT VIEPSQKAMK ESVRFLTDFG EISDAISDLL TSSPNFNVIS AIGPQGAGKS TLLSMLAGNN SRQMYREYVF RVSREANEQ SRHQTIQIDI YIVNHQIFLD CQPMYSFSIM EGLPKVRGGR FDDSTAMSDT LRLTAFLLYV SHTVLVVSET HYDKVIIDTL RVAEQIRPYL AIFRPKLAID RKTNLVFIKT KASSIDLAPT VIREREELLR LSFQDSRWLK VSQEPFKTLI VLEEIRVRRE HLFEEGDEPD EAASLNEFDE QIAELREELQ KNREDFTVET AAMDEKKWLD MCREVIRDKT LHKTLLKEYQR AMTDGVRTHF DNGFHAERDA NKFFS
Source	E.coli
Target Names	smg-9
Protein Names	Recommended name: Protein smg-9 Alternative name(s): Suppressor with morphogenetic effect on genitalia protein 9
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.