



Recombinant *Drosophila melanogaster* Protein cutoff (cuff)

Product Code	CSB-MP896050DLU
Abbreviation	cuff
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V629
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNSNYTILNI RAHSWAEKDN LVFPTITKPV SNGWFAWSPI GQFEPNSQKV PYVVPPSIK FPISLRYKDS PPKLEKKPNI FLDNMLQYIE SSSYMYVMVR NNQQAQLNAN IVCTSEVLEL MMCAPYEKKT GWSLGVTRYR NTMYICRIDV EQDPIDQDN LKRAMQEFWL RNLRTHCVFE NGIKMHQHNQ SSEEHLRFHG VFSFDLNGNR VLFDSPVLAE MPSTTLNGSA LSWVDLQIRP MFMSRLDWPE HNRTEALKWW VKCFLLGIES LYIARRDENA HVHNIKTLV RDLWKSCEKD WSPTVCANFM IYLLNCISQV MAPIDCPSTV YLFQFDASQG TVSYKGLRGR NQYTFVSDWF RMMLDDHTND MCKTPNLQTM SSIV
Source	Mammalian cell
Target Names	cuff
Protein Names	Recommended name: Protein cutoff
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.