



Recombinant Drosophila melanogaster SH3 domain-binding protein 5 homolog (pcs)

Product Code	CSB-MP896054DLU
Abbreviation	pcs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V785
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSSAEDGELD PQIQIELENL NSATDEINKL EIELEEANST FRILLNESTR RLKVSSKKLG NCIEKARPYE EALDKAREAQ IECQKAAVKF QRANEIHAAA KETVALAEQR FMSNSHEWQF DNAWQEMLNH ATQKVMDAET QKADCHAEHQ RLTKLFNAAE QKLQQLEDRF RRSINKSRPY FEEKQVCQDQ LQTQKNRIQE LQQQVAGAKS TYSTALRNLE RISEDHRQR GDFPTPPGPR EPGVGAELNS PTSSALPSLP DFQLELEKCD YPSIAGSQMS LGAKTPQAAA ETEDEEDACD YDETGAGELR GVVDERDLEA LRQVKILAV RPIEGGDGQQ QNDVWEHELK ATVDKLDHLM MLKETAKRQQ TNRLKSTEQR PDSLGAEALK RHCDVVEVKV TSCATTASLP VTPHHQLNHL APPTPIKKLQ QQLAPLPSVN VSMRELPLLA RLSNELLD RS SAAF GGVRKT LRRRSLE
Source	Mammalian cell
Target Names	pcs
Protein Names	Recommended name: SH3 domain-binding protein 5 homolog Alternative name(s): Protein parcase
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.