



Recombinant Drosophila melanogaster S-methyl-5'-thioadenosine phosphorylase (CG4802)

Product Code	CSB-MP896056DLU
Abbreviation	CG4802
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V813
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MITIKCKDTN LDPLPVKIGI IGGSGLDLDDPD ILEQRQERVV ETPYGEPSDA LIEGEINGVQ CVLLARHGRK HDIMPNSVNY RANIWALRDV GCTHLIVSTA CGSLREEIKP GNLVVPDFI DRITKRLQTF YDGKAQSPRG VCHLPMFPAF SERTRNILLQ AAKELEIPAH DKATIVTIEG PRFSSRSESH MFRQWGGDLI NMTTCPEVVL AKEAGLLYGS VAIATDYDCW RMGCEGVNVQ DVLRTFAENV IKVKILVNA VGRIAKEDWS EDILNAKQCV CNNTMSGAM
Source	Mammalian cell
Target Names	Mtap
Protein Names	Recommended name: S-methyl-5'-thioadenosine phosphorylase EC= 2.4.2.28 Alternative name(s): 5'-methylthioadenosine phosphorylase Short name= MTA phosphorylase Short name= MTAP Short name= MTAPase
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.