



Recombinant *Drosophila melanogaster* Methylthioribose-1-phosphate isomerase (CG11334)

Product Code	CSB-EP896061DLU-B
Abbreviation	CG11334
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V9X4
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSLQSIKYSR GSLEILDQLL LPGQSKYVVV RGVEDGWKVI NKMQVRGAPA IAIVGCLSLA VEINPEDFEN KKSLRQEIEG KLNLYLSARP TAVNMKIAAD ELITLANELY KDEAIDVTQM KHRFLDATEA MLKKDIADNR AIGANGAQAI LQGVAKAGKA PAGSTGSVRV LTHCNTGSLA TAGYGTALGV VRQLAELGKL EHVYCTETRP YNQGARLTAY ELVHEKFPAT LVLDSMVAAL LRAKNVAADV VGADRVASNG DTANKIGTYQ IAVVAKHHDV PFYVAAPLTS IDLAIPGGDH IIIEERPDRE MTHVGEHRIA APGINCWNPV FDVTPASLIT GIITERGVFK PAELKEAITK LLES
Source	<i>E. coli</i>
Target Names	CG11334
Protein Names	Recommended name: Methylthioribose-1-phosphate isomerase Short name= M1Pi Short name= MTR-1-P isomerase EC= 5.3.1.23 Alternative name(s): S- methyl-5-thioribose-1-phosphate isomerase Translation initiation factor eIF-2B sub
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.