



Recombinant Drosophila melanogaster Probable DNA-directed RNA polymerase III subunit RPC6 (CG5380)

Product Code	CSB-EP896075DLU
Abbreviation	CG5380
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VD25
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MATEVSQLLL AVVQGIPAGA TNDDLTKALP DVPAATRVEA LNILLQEGGI EILKKGEKLV YRAKDPEKKS ALPKDADNEE KVVYGIVEEG GNKGIWIRDI RMKSNLNIQ LNKILKNLET KKLKAVKSV NASKKKVYML YNLEPDLISIT GGAWYQDQDF EVEFVDVLNQ QCLRFLQMKR DSAEKKREGP LAFKQMSCCT VNEVQKFISD LGISKVNLAE ADLETILKTV VYDGNAERVR QQDGSFVYRA VNAPLPPTGL VQMPCGICPV IKNCSNCGDV TAITCEYMRD WLD
Source	E.coli
Target Names	CG5380
Protein Names	Recommended name: Probable DNA-directed RNA polymerase III subunit RPC6 Short name= DNA-directed RNA polymerase III subunit F Short name= RNA polymerase III subunit C6
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.