



Recombinant *Drosophila melanogaster* Protein BCCIP homolog (CG9286)

Product Code	CSB-BP896088DLU
Abbreviation	CG9286
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VFR0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSANKQKKLS TMEVDPNEDV SSSSEDDDDD EPHPDAYKGN EEVQIDFEGR APVDPDAQGI SLLQLRFLR AHINCNQMAD LIIAQNFIGS VICQCDDEGT ESETEDDNMV EDGTIFGITS VLNLTAKKDQ PSIAQLRTYI LDRAKTHASP EVQQQLKEIL DSEQRHVGFL INERFINIPA QISVPLLQSL QQEIEAAKAK KMKYDFGTLL LLVKFYRKES KKGKPGEDVY TNAEDELSD RAKFSFEYSM ASETDSGMSG DWLEGDAVMT PYRKLVLVLEA KKLPQLIDDI QRFINGE
Source	Baculovirus
Target Names	CG9286
Protein Names	Recommended name: Protein BCCIP homolog
Expression Region	1-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.