



# Recombinant Drosophila melanogaster Zinc carboxypeptidase A 1 (CG17633)

<b>Product Code</b>	CSB-EP896112DLU-B
<b>Abbreviation</b>	CG17633
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9VL86
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila melanogaster (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ERVRYDNY RMYKVNSENA KQLEVLKDLE GSSDSIMFLD GVHLVGADIQ IIVAPHKVPD FLEILGKSEI KYELQSRDVQ KSLDEIDEKV AIKGRATTAY NWAQYYELDD TYAWLQSLAQ TNPGVVTLIE GGKTYQGRSI LGVKITKGGE TINGKAKPGI FLEAGIHARE WIAPAAATFI INQLLTSEVE NIKELAENYT WYVLPHPANPD GYVYTHHTNR LWRKTRTPYG SCFGADPNRN WGFHWNEVGA SSSACSDTYA GPSAFSEIET LSLSKFIEGL KGKVQLYLSL HAYSQYLLYP YGHTSDLPDN VADFEKVFDA SIAAVNKRYG TTYTGGNIYD AIYPAAGASV DWAYGTQDVR MAFCYELRPS STSYLTGFKL PAEQIVPASE ELLDSIVAMA TEVKSLGYFD
<b>Source</b>	E.coli
<b>Target Names</b>	CG17633
<b>Protein Names</b>	Recommended name: Zinc carboxypeptidase A 1 EC= 3.4.17.-
<b>Expression Region</b>	23-430
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.