



Recombinant *Drosophila melanogaster* Uncharacterized protein CG9548 (CG9548)

Product Code	CSB-EP896119DLU-B
Abbreviation	CG9548
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VMC8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAKHHPDLIF CRKQPGVAIG RLCEKDDGKC VICDSYVRPC TLVRICDECN YGSYQGRCVI CGGPGVSDAY YCKSCTIQEK DRDGCPKIVN LGSSKTDLFY ERKKYGFKQN Y
Source	<i>E. coli</i>
Target Names	Phf5a
Protein Names	Recommended name: Uncharacterized protein CG9548
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.