



Recombinant Drosophila melanogaster Protein extra bases (exba)

Product Code	CSB-EP896125DLU-B
Abbreviation	exba
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VNE2
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSQKTERPVL SGQRIKTRKR DEREKYDPTG FRDAVIAGLE KTEGDLEQIS KYLDSAGNKL DYRRYGEVLF DILIAGLLV PGGISISQDGE KPRTSYCIFD APESMESMRN HEQVFVKLIR RYKYLEKMFE EEMGKVLLFV KGFTPSEIRIK LARMTALWLV NGSVPPNVLL VLNNEHLIKD GIALEFLEL FQTFKQEKGI AYLIQALKKG GLESKLMDFF PPNKRTEEYF KQVFLDKELN EIVKLHKAQA SQEAKRELQQ ALIDDINDEK PYNEITSDIK DFSQRTNIPD HEIIVIIWST IMSLGEWNKK EELVTDQAVR HLKNYCPLLQ AFASTDRSEL ALILKVQEFC YENMNFMKAF QKIILLFYKT EVLSEEIILR WYKEGHSNKG KMHFLEQMRK FVEWLQSAEE ESESEDEQKN GE
Source	E.coli
Target Names	kra
Protein Names	Recommended name: Protein extra bases Alternative name(s): Eukaryotic translation initiation factor 5C domain-containing protein Short name= Decp Short name= eIF5C domain-containing protein
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.