



Recombinant Drosophila melanogaster Putative N-acetylglucosamine-6-phosphate deacetylase (CG17065)

Product Code	CSB-EP896136DLU
Abbreviation	CG17065
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VR81
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDVRIKCNSE EPESPEQKIL ASSQRLQFT NCRLVRDHRI IHDDLWVRDG RIVNPEPVFF DERTKAHCRI DCGGAIAPG YIDLQINGGY GVDFSYDTET IEEGVATVAR GLVKSGVTSF CPTLVTSFSD SYHTILPRIP AEVPKGAGIL GIHAEGPFIN PQKKGAPHEH CIQTIDKGLS TLKETYGSLE RIKIITLAPE KVTDPVEVIGQ LVERGITVAL GHSMASLSDG ERAVQQGATL ITHLFNAMLP FHHRDPGLVG LLASDAVPHG RTVYFGIISD GVHTHPAALR IAYRTHPQGL ILVTDAISAL GLEEGVHHIG QLPLQVKQKG AFIAGTETLC GSIAPMDECV RIFQKATDCS VVYAIEAATL HPAQCLKIEK QKGTLDGSD ADFVLLDDQL RVLSTWIAGV CVHRTVK
Source	E.coli
Target Names	CG17065
Protein Names	Recommended name: Putative N-acetylglucosamine-6-phosphate deacetylase EC= 3.5.1.25 Alternative name(s): GlcNAc 6-P deacetylase
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.