



Recombinant Drosophila melanogaster Nitric oxide synthase-interacting protein homolog (CG6179)

Product Code	CSB-EP896161DLU-B
Abbreviation	CG6179
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9VWV8
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MTRHARNCTA GAVYTYNEKK RDAAESGYGT NAQRLGKDSV KSFDCCSLTL QPCRRPVITK DGYLFDKEAI LQYIVTKKNE YSRRLKEYER LRRAEEDKLS QEANSKQQR MERFVNAEKP AMTPAHSSAA ASEKPSTSSA AAAASSESSS ASSISNMTNG HEKKLPSFWL PSECPNAGLA KAQKPDATYI CPVSQKPLRV KDLIDVKFTL LKDGDTKRSL IAKEARYMCP ITHDVLSNAV PCAVLRPTGD VVTMECVERL IRKDMIHPLT DRKLKEKDII PLQRGGTGYA TTNDHLQAKE KRPMLQA
Source	E.coli
Target Names	CG6179
Protein Names	Recommended name: Nitric oxide synthase-interacting protein homolog
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.