



Recombinant *Drosophila melanogaster* Mediator of RNA polymerase II transcription subunit 30 (MED30)

Product Code	CSB-MP896176DLU
Abbreviation	MED30
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9W0P3
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MWKYGQNQGN QGPSSGGGGG GGPNMMPMGG FGMQHGNMQQ MHMSPQHQQQ QQQMGMGGP GSMQMNPQGP GPGGLMPGM SPQHMQQQQ QQQMMQQMM VPQQGVGVGM GGGGVGMGG GGVVPQQQQ QPQQNMPQQN IPQQQQQLNP VAGIPPGGAG GSNNMLAISQ QNPHEINIV QLSRLGQETV QDIASRFQEV FASLKGIQPT SHRENSSEKK VQEYFRTIRL LFKRVRIIYE KCNDAGMDYM SAESLIPYRD EPEPRIEPSL CDEYRKVLQE NHELIVTVKL KNRQLREIID RTRIIIWEIN TMLAMRRS
Source	Mammalian cell
Target Names	MED30
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 30 Alternative name(s): Mediator complex subunit 30 dMED20 dTRAP25
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.