



Recombinant Human immunodeficiency virus type 1 group M subtype J Gag polyprotein (gag)

Product Code	CSB-BP896215HHAD
Abbreviation	gag
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WC53
Product Type	Recombinant Protein
Immunogen Species	Human immunodeficiency virus type 1 group M subtype J (isolate SE9280) (HIV-1)
Purity	≥85% (SDS-PAGE)
Sequence	GARASILSG GKLDWEKIR LRPGGKKKYR IKHLVWASRE LDRFALNPGL LESAGCQQI LVQLQPALQT GTQEIKSLYN TVATLYCVHQ RIEIKDTMEA LEKIEEIQNK NKQQAQKAET DKKDNSQVSQ NYPIVQNLQG QPVHQALSPR TLNAWVKVIE EKAFSPEVIP MFSALSEGAT PQDLNTMLNT IGGHQAAMQM LKDTINEEAA EWDRVHPVHA GPIAPGQVRE PRGSDIAGTT STLQEQIGWM TGNPPIPVGE IYKRWILGL NKIVRMYSVP SILDIRQGPK EPFRDYVDRF FKALRAEQAT QDVKNWMTDT LLVQANPDC KTILKALGSG ATLEEMMTAC QGVGGPGHKA RVLAEAMSQV TNTNIMMQRG NFRDHKRIVK CFNCGKQGGHI AKNCRAPRKK GCWKCCKEGH QMKDCTERQA NFLGKIWPSS KGRPGNFLQS RPEPTAPPAE SLGLGEEIPS PKQEPKDKEL YPLTSLKSLF GSDPLSQ
Source	Baculovirus
Target Names	gag
Protein Names	Recommended name: Gag polyprotein Alternative name(s): Pr55Gag Cleaved into the following 6 chains: 1. Matrix protein p17 Short name= 2. MA 3. Capsid protein p24 Short name= 4. CA 5. Spacer peptide p2 6. Nucleocapsid
Expression Region	2-497
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.